<image>

COYOTE[®] FlashDetect™ LyocartE MRSA/SA assay

Effective differentiation of MRSA and SA in less than 35 mins

All-in-One Cartridge

Less hands-on time, less operational error

FlashDetect™

- Differentiate between MRSA and SA in one reaction
- Sample to result within 35 mins
- LoD 500 copies/mL
- Includes internal control
- Pre-loaded with lyophilized reagents
- Contamination-free
- Room temperature storage and transportation

Compatible with Multiple Clinical Specimens

- Nasal swabs
- Skin and soft tissue swabs
- Blood culture samples






TECHNICAL SPECIFICATIONS

FlashDetect™ LyocartE MRSA/SA Assay

····· · · · · · · · · · · · · · · · ·				
Catalog #	6018020305-CE	6018020305-CE		
Method	Real-Time PCR			
Fluorescence Channel	FAM	ROX		
	Staphylococcal protein A (SPA) gene	Staphylococcal cassette chromosome mec (SCCmec) gene		
	HEX	CY5		
	mecA/mecC gene	RNase P gene		
LoD	500 copies/mL	500 copies/mL		
Detection Time	35 mins	35 mins		
Hands-on Time	< 1 min	< 1 min		
Sample Type	Nasal swabs, skin and soft tis	Nasal swabs, skin and soft tissue swabs, blood culture samples		
Integrated workflow	, , , , , , , , , , , , , , , , , , ,	Fully automated process for extraction, amplification and result interpretation with FlashDetect™ system		
Compatible Instrument	FlashDetect™ Flash10 and Na	FlashDetect™ Flash10 and Nano		
Storage Temperature	2-28°C			



WORKFLOW



ORDERING INFORMATION

Compatible Instrument

Model	Throughput	Dimensions (L x W x H)
FlashDetect™ Nano	1	12.5 x 22 x 29.5 cm
FlashDetect™ Flash10	4	38 x 28.7 x 35 cm

Cartridge

Product	Pack Size	Catalog #
FlashDetect™ LyocartE MRSA/SA assay	24 tests	6018020305-CE



Coyote Bioscience Co., Ltd.

H.O.: Building 22, Zone 3, Gaolizhang Road, Haidian District, Beijing 100095, China. Tel: +86 10 64844237 Email: sales@coyotebio.com