

VIRPLEX MiniLab

SYNDROMIC POINT-OF-CARE MOLECULAR TESTING

- RealTime PCR system with integrated extraction.
- Simultaneous detection capacity of up to 24 targets.
- Closed-loop process with no risk of contamination.
- From sample to result in less than 2 h.
- Minimum hands-on time (< 2 minutes).
- User-friendly graphical interface with Ct value and automatic result analysis.
- Enhance your laboratory performance in the most space-efficient way.



DESCRIPTION	CAT. No.	PACK SIZE
MENINGITIS VIRPLEX MINILAB PANEL*	ML002	10 test
STD VIRPLEX MINILAB PANEL* <i>Chlamydia trachomatis</i> , <i>Chlamydia trachomatis</i> LGV, <i>Neisseria gonorrhoeae</i> , <i>Trichomonas vaginalis</i> , <i>Mycoplasma genitalium</i> (MGE), MGE macrolide resistance, <i>Ureaplasma urealyticum</i> , <i>Ureaplasma parvum</i> , <i>Mycoplasma hominis</i> , Herpes simplex virus 1, Herpes simplex virus 2, <i>Haemophilus ducreyi</i> , <i>Treponema pallidum</i> , Monkeypox virus, Internal control	ML001	10 test

*Available soon